**BD**   **Phoenix™**  
**NMIC/ID Panel**

**INTENDED USE**

The **Phoenix™** NMIC/ID panel is used for the rapid identification and susceptibility testing of most aerobic and facultative anaerobic gram-negative bacteria of human origin with selected antimicrobial agents (for a complete listing of taxa, refer to the **Phoenix** System User's Manual). This panel is only for use with the **Phoenix** Automated Microbiology System instrument (Cat. No. 448100).

**PRINCIPLES OF THE PROCEDURE**

Refer to the **Phoenix** System User's Manual.

**PANELS**

*Identification*: Refer to the **Phoenix** System User's Manual.

*Susceptibility*: The NMIC/ID panel contains the antimicrobial agents and concentrations in doubling dilutions found in Table 1.

**STORAGE**

Store at 15–25° C. Do not use panel if the pouch is punctured or opened. Do not use the panel if desiccant is missing or if the desiccant pouch is torn.

Panels must be used within 2 h of being removed from the pouch.

**Precautions**: For *in vitro* Diagnostic Use.

**WARNINGS**

Observe established precautions against microbiological hazards throughout all procedures. "Standard Precautions"<sup>1,2</sup> and institutional guidelines should be followed in handling all items contaminated with specimens and microorganisms. Prior to discarding, sterilize specimen containers and other inoculated materials by autoclaving.

**QUALITY CONTROL**

*Identification*: See expected results below.

Test Organism	Expected Result
<i>Escherichia coli</i> ATCC™ 25922	<i>Escherichia coli</i>
<i>Pseudomonas aeruginosa</i> ATCC 27853	<i>Pseudomonas aeruginosa</i>

*Susceptibility*: See expected results in Table 1.

NOTE: QC may be reported as less than or equal to the lowest, or greater than the highest concentration of the antimicrobic.

**LIMITATIONS OF THE PROCEDURE**

The clinical relevance of a specific antimicrobial and organism combination is determined by the activity of that antimicrobic against the organism and whether the antimicrobic is indicated for treatment of a disease state associated with that organism. The **Phoenix** System provides results for combinations, whether they are clinically relevant or not relevant. Overall, the results obtained using the **Phoenix** System compare favorably with the NCCLS broth microdilution reference methodology<sup>3</sup>; however, Table 2 contains clinically relevant combinations that did not meet the strict standards of BD during clinical trials and will either not be reported or an alternate method is recommended for confirmation of the result. Similarly, Table 3 contains combinations which are not clinically relevant and are either not reported or an alternate method is recommended for confirmation of the result.

**PERFORMANCE CHARACTERISTICS**

Refer to the **Phoenix** System User's Manual.

**REFERENCES**

- National Committee for Clinical Laboratory Standards. 2001. Approved guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2<sup>nd</sup> ed. NCCLS, Wayne, Pa
- U.S. Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030. 1991. Occupational exposure to bloodborne pathogens. Federal Register 56:64175-64182.
- National Committee for Clinical Laboratory Standards. 2003. Approved standard M7-A6. Methods for dilution antimicrobial susceptibility tests for bacteria that grow aerobically, 6<sup>th</sup> ed. NCCLS, Wayne, Pa.

ATCC is a trademark of the American Type Culture Collection. BD, BD Logo and Phoenix are trademarks of Becton, Dickinson and Company. ©2004 BD.

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